



Affordable Reliability

Connect the Big-Lifter's® docking device to a scooter or powerchair, push a button on the hand-held pendant, and the Big-Lifter will automatically raise the mobility device to the desired height. Manually rotate the scooter or powerchair into the vehicle, press a button and the Big-Lifter will gently lower the mobility device into the rear cargo area.

Fully Powered. Automatic lifting and lowering of a mobility device.

Compact Design. Partial third-row seating or extra cargo space is available when a mobility device is not present (optional quick-release pins required).

Cost-Effective. Manual rotation of the mobility device enhances the Big-Lifter's affordability.

Safe. Meets or exceeds all Federal Motor Vehicle Safety Standards.

Adaptable. Reinstall on another applicable vehicle in the future.

Lifting Capacity

400 lb (181 kg)

Vehicles


SUVs, Vans

Warranty

3-Year Limited

Learn More

Bruno.com



OPTIONS

- Back-Off Device (A)
Removes back of mobility device to fit in small openings
- Driver Side Placement (B)
- Exterior Hand-Held Control Package (C)
- Quick-Release Pin (D)
Allows lift head to fold down

